


Testing For Endotoxins?

Associates of Cape Cod, Inc. offers a USP BET¹ compliant endotoxin testing system!



For fully USP BET¹ compliant endotoxin tests, use ACC's new Pyros Kinetix[®] Flex

USP Bacterial Endotoxins Test¹ Requirements for Photometric Quantitative Techniques (Test for Interfering Factors/Test Procedure)

USP BET Requirements	Pyros Kinetix [®] Flex Tube Readers
A. Sample solution (without added endotoxin), ≥ 2 replicates	✓
B. Sample solution with endotoxin added to give a concentration from the middle of the standard curve, ≥ 2 replicates	✓
C. At least three standard endotoxin concentrations, ≥ 2 replicates of each	✓
D. Water for BET (negative controls), ≥ 2 replicates	✓

Reference:

¹Bacterial Endotoxins Test. Chapter <85>. USP 35, p. 88-92, 2011, United States Pharmacopeial Convention, Rockville, MD, USA.



ASSOCIATES OF
CAPE COD
INCORPORATED

Specialists in Endotoxin and Glucan Detection
888.395.2221 • WWW.ACCIUSA.COM



Incubating Kinetic Tube Reader



ISO 13485:2003
MEDICAL DEVICES