


# Testing For Endotoxins?

**Associates of Cape Cod, Inc. offers a USP BET<sup>1</sup> compliant endotoxin testing system!**



## For fully USP BET<sup>1</sup> compliant endotoxin tests, use ACC's new **Pyros Kinetix<sup>®</sup> Flex**

### USP Bacterial Endotoxins Test<sup>1</sup> Requirements for Photometric Quantitative Techniques (Test for Interfering Factors/Test Procedure)

USP BET Requirements	Pyros Kinetix <sup>®</sup> Flex Tube Readers
A. Sample solution (without added endotoxin), $\geq 2$ replicates	✓
B. Sample solution with endotoxin added to give a concentration from the middle of the standard curve, $\geq 2$ replicates	✓
C. At least three standard endotoxin concentrations, $\geq 2$ replicates of each	✓
D. Water for BET (negative controls), $\geq 2$ replicates	✓

Reference:

<sup>1</sup>Bacterial Endotoxins Test. Chapter <85>. USP 35, p. 88-92, 2011, United States Pharmacopeial Convention, Rockville, MD, USA.