

Fagron SyrSpend® SF



Superior Suspending Technology

Compounding pharmacists need to have quality stability and BUD data to support their efforts creating customized medications. Over 70 APIs with SyrSpend® SF products have been studied and published in reputable peer review journals, including IJPC. Five of the top 10 drugs administered to 38 million beneficiaries enrolled in Medicare Part D have stability and BUD supporting data.¹ To learn more about the APIs studied with SyrSpend® SF products, visit us.fagron.com/en-us/syrSpend.

SyrSpend® SF Vehicles

- Extensive and growing database of scientific references for stability data and beyond-use dating
- Pre-weighed convenience packs available to compound suspensions in minutes
- Supports adjusted and flexible dosage regimens

SyrSpend® SF

- All-in-one syrup suspending vehicle
- Simple one-step compounding process
- pH adjusted to 4.2 for maximum compatibility
- Also available as powder for reconstitution (preservative-free)

SyrSpend® SF Alka

- Easy-to-use suspending powder for reconstitution
- Provides an alkaline (>7) environment for acid-labile drugs
- Preservative-free

Active Suspending Technology

- Holds particles in suspension longer
- Allows more accurate dosing

Low Osmolality

- <50 mOsm/kg
- Minimizes GI upset and diarrhea
- No sorbitol means lower osmolality to lessen laxative/diarrhea problems
- Ideal for pediatric and geriatric patients

Sweetened without Sugar or Alcohol

- Bitterness is masked
- Pleasant taste enhances compliance



¹<https://www.cms.gov/Newsroom/MediaReleaseDatabase/Fact-sheets/2016-Fact-sheets-items/2016-08-18.html>

fagron.us

Fagron, Inc. - 2400 Pilot Knob Road - St. Paul, Minnesota 55120 - USA
Tel. (800) 423 6967 - Fax (800) 339 1596 - info@fagron.us

