

WAY

14

MEDISCA Presents the MD3™ Multiple Dosages. One Ground-Breaking Solution.

The MD3 is an innovative multidose actuator that delivers numerous benefits to the well-known triad relationship between patients, pharmacists, and physicians. With three doses in one, the MD3 represents the next step in medical dispensers.

President & CEO: Antonio Dos Santos

Founded: 1989

Employees: 400

Toll-free Phone: (800) 932-1039

Fax: (855) 850-5855

Address: MEDISCA 661 Route 3, Unit C
Plattsburgh, NY 12901, USA

Website: www.medisca.com



Company Background

With over 25 years of experience, MEDISCA has developed a solid reputation as an innovative and entrepreneurial company in the pharmaceutical compounding industry. MEDISCA is a leading FDA-registered supplier of quality pharmacy compounding products and offers excellent service to the North American, Australian, and international markets. MEDISCA was recently recognized with an ISO 9001:2008 certification from TUV, demonstrating our strong commitment to the quality of our products, the satisfaction of our customers, and the continuous improvement of our internal processes.

Product Overview

The MD3 is a revolutionary patent pending actuator that can deliver three different dosages while providing extreme precision and accuracy. This revolutionary device is the newest addition to the MD line of dispensing products and delivers numerous benefits to the well-known triad relationship between patients, pharmacists, and physicians.

Features & Options

With three doses in one, the versatile delivery system of the MD3 effectively replaces the need for three separate dispensers. This revolutionary actuator provides a number of great features and benefits, such as:

- Multiple dosage options (0.1, 0.25, and 0.5 mL).
- Patented airless dispensing technology.
- Ideal for fluctuating prescription regimes.
- Compatible with MEDISCA's entire line of UV-resistant MD bottles ranging from 30 mL to 240 mL.
- Elegant design with multiple color options such as green, blue, and pink.
- Accuracy and precision results available.

Testimonial

“With multiple dispensing options, the MD3 has revolutionized the way dosages can be managed. This product is ideal for clinical applications that require dosages to be titrated. This will optimize symptom control for a variety of clinical conditions such as hormone restoration therapy and transdermal pain treatments. Patients and caregivers will benefit from the simple design, physicians can better manage medication regimens, and pharmacists like myself will no longer require multiple dispensers. Its ease of use promotes patient compliance.”

— Ken Speidel, RPh., Pharm D.

Three Simple Steps for Dispensing



1. SET THE DOSAGE

Turn the actuator clockwise to desired dosage.



2. DISPENSE

Immediately turn counter clockwise. Return to its original position.



3. APPLY

Apply the medication. No need for application pad.

Visit MD3.com to watch a video highlighting the many features and benefits of the MD3.