



Superior Accuracy for Superior Results



WATCH OUR NEW VIDEO!



Snapping made
quick and effortless!

Revolutionizing Airless Dispensing

MD CapSnapper

New Items Based On Popular Demand

PRODUCT	SIZE	DOSAGE	ACCURACY (%)
 MD Pen	6.5mL	 0.10mL, 0.15mL	+/-2%
 MD Syringe	5mL, 10mL	0.15mL	+/-2%
 MD Jar	50mL	0.5mL	+/-10%
 MD Pump	30mL, 50mL, 80mL, 100mL, 150mL, 200mL,  240mL	Actuator A: 0.5mL Actuator B: 1.0mL  Actuator C: 1.5mL	+/-5%

ORDER NOW!



1.800.932.1039

| WWW.MEDISCA.COM

